


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 4.9mm

180 gms: angled; fenestrated Ø 3.5mm



SKU: 244-433-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.58

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.58

SPECIFICATIONS

[angled](#), [permanent](#)