


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 9mm; opening 9.1mm

150 gms; straight; fenestrated Ø 5mm




**SKU:** 244-436-09T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.59.61
<b>SPECIFICATIONS</b>	<a href="#">permanent</a> , <a href="#">straight</a>

**Schließkraft**  
**Closing Force**



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.59.61

## SPECIFICATIONS

[permanent](#), [straight](#)