


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 12mm; opening 10.3mm

180 gms; straight; fenestrated Ø 5mm




**SKU:** 244-436-12T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                 |
| <b>TYPE</b>           | <a href="#">Vessel Clips</a>                         |
| <b>MATERIAL</b>       | <a href="#">Titanium</a>                             |
| <b>AUTHORS</b>        | <a href="#">YASARGIL</a>                             |
| <b>Medicon</b>        | 58.59.62   |
| <b>SPECIFICATIONS</b> | <a href="#">permanent</a> , <a href="#">straight</a> |

**Schließkraft**  
**Closing Force**



### ZUSÄTZLICHE INFORMATIONEN

|                 |                              |
|-----------------|------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>         |
| <b>TYPE</b>     | <a href="#">Vessel Clips</a> |
| <b>MATERIAL</b> | <a href="#">Titanium</a>     |
| <b>AUTHORS</b>  | <a href="#">YASARGIL</a>     |
| <b>Medicon</b>  | 58.59.62                     |

## SPECIFICATIONS

[permanent](#), [straight](#)