


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 5mm; opening 5.5mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-439-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.66
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.66

SPECIFICATIONS

[angled](#), [permanent](#)