


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 11mm; opening 5.5mm

180 gms: angled; fenestrated Ø 5mm



**SKU:** 244-439-11T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.59.68
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.59.68

## SPECIFICATIONS

[angled](#), [permanent](#)