


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 20mm; opening 11.4mm

200 gms; straight




SKU: 244-400-20T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.21
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.21

SPECIFICATIONS

[permanent](#), [straight](#)