


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 5.8mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-443-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.78

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.78

SPECIFICATIONS

[angled](#), [permanent](#)