


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 7.8mm

150 gms; straight; fenestrated Ø 5mm




SKU: 244-445-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.59.82 |
| SPECIFICATIONS | permanent , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.59.82 |

SPECIFICATIONS

[permanent](#), [straight](#)