


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 5mm; opening 5.8mm

150 gms; angled; fenestrated Ø 5mm




SKU: 244-449-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	angled , permanent

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	angled , permanent