


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.2mm

150 gms; angled; fenestrated Ø 5mm




SKU: 244-449-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.89
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.59.89

SPECIFICATIONS

[angled](#), [permanent](#)