


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6mm

150 gms; curved



SKU: 244-401-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)


Aesculap FT742T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.08

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT742T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.61.08

SPECIFICATIONS

[curved](#), [permanent](#)