

## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 7mm

90 gms; straight



**SKU:** 244-500-09T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT250T

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.60.60

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT250T

**MATERIAL** [Titanium](#)

## AUTHORS

[YASARGIL](#)

## Medicon

58.60.60

## SPECIFICATIONS

[straight](#), [temporary](#)