


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 11mm; opening 7.8mm

90 gms; straight



SKU: 244-500-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT100T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.62
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT100T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.60.62

SPECIFICATIONS

[straight](#), [temporary](#)