


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 11mm; opening 7.8mm

90 gms; straight



**SKU:** 244-500-11T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT100T
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.60.62
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">temporary</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FT100T
<b>MATERIAL</b>	<a href="#">Titanium</a>

## AUTHORS

[YASARGIL](#)

## Medicon

58.60.62

## SPECIFICATIONS

[straight](#), [temporary](#)