


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 15mm; opening 9.2mm

90 gms; straight



SKU: 244-500-15T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|------------------------------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.60.66 |
| SPECIFICATIONS | straight , temporary |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.60.66 |

SPECIFICATIONS

[straight](#), [temporary](#)