


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 25mm; opening 13.3mm

110 gms; straight



SKU: 244-500-25T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.60.72 |
| SPECIFICATIONS | straight , temporary |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.60.72 |

SPECIFICATIONS

[straight](#), [temporary](#)