


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 25mm; opening 13.3mm

110 gms; straight




SKU: 244-500-25T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.72
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.60.72

SPECIFICATIONS

[straight](#), [temporary](#)