


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 6mm

110 gms; curved



SKU: 244-501-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.61.58 |
| SPECIFICATIONS | curved , temporary |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.61.58 |

SPECIFICATIONS

[curved](#), [temporary](#)