


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 20mm; opening 11mm

90 gms; curved




SKU: 244-501-20T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.71
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.71

SPECIFICATIONS

[curved](#), [temporary](#)