


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 5.5mm

110 gms; curved



SKU: 244-503-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.72
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.72

SPECIFICATIONS

[curved](#), [temporary](#)