


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 15mm; opening 7.4mm

110 gms; curved



SKU: 244-503-15T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.86
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.86

SPECIFICATIONS

[curved](#), [temporary](#)