


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.2mm

151 gms; straight



**SKU:** 244-400-07T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT740T

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.60.08

**SPECIFICATIONS** [permanent](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT740T

**MATERIAL** [Titanium](#)

## AUTHORS

[YASARGIL](#)

## Medicon

58.60.08

## SPECIFICATIONS

[permanent](#), [straight](#)