


## YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 6.2mm

150 gms; straight



**SKU:** 244-400-07

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE740

**AUTHORS** [YASARGIL](#)

**SPECIFICATIONS** [permanent](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE740

**AUTHORS** [YASARGIL](#)

**SPECIFICATIONS** [permanent](#), [straight](#)