

ZENKER RETRACTOR

75x12mm;; 35° angled

27.5cm (11")


SKU: 064-320-27

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	BT596
KLS Martin	15-321-03
AUTHORS	ZENKER
Medicon	23.29.75
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	BT596
KLS Martin	15-321-03

AUTHORS

[ZENKER](#)

Medicon

23.29.75

SPECIFICATIONS

[angled](#)