

## ZENKER RETRACTOR

75x12mm;; 35° angled

27.5cm (11")

**SKU:** 064-320-27

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Retractors</a>
<b>Aesculap</b>	BT596
<b>KLS Martin</b>	15-321-03
<b>AUTHORS</b>	<a href="#">ZENKER</a>
<b>Medicon</b>	23.29.75
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Retractors</a>
<b>Aesculap</b>	BT596
<b>KLS Martin</b>	15-321-03

## AUTHORS

[ZENKER](#)

## Medicon

23.29.75

## SPECIFICATIONS

[angled](#)