

ZENKER RETRACTOR

35 x 8mm; 35° angled

25.0cm (10")


SKU: 064-320-25

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	BT594
KLS Martin	15-321-01
AUTHORS	ZENKER
Medicon	23.29.35
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	BT594
KLS Martin	15-321-01

AUTHORS

[ZENKER](#)

Medicon

23.29.35

SPECIFICATIONS

[angled](#)