

## ZENKER RETRACTOR

35 x 8mm; 35° angled

25.0cm (10")


**SKU:** 064-320-25

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Retractors</a>
<b>Aesculap</b>	BT594
<b>KLS Martin</b>	15-321-01
<b>AUTHORS</b>	<a href="#">ZENKER</a>
<b>Medicon</b>	23.29.35
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Retractors</a>
<b>Aesculap</b>	BT594
<b>KLS Martin</b>	15-321-01

## AUTHORS

[ZENKER](#)

## Medicon

23.29.35

## SPECIFICATIONS

[angled](#)